
# MultiMode AFM Fluid Precautions

## 246.1 Introduction

This support note describes a service procedure for use with Digital Instruments SPM's performing fluid imaging. Users who are performing fluid imaging should conduct the inspection and maintenance routine described below on a regular basis. This, combined with spillage precautions, will extend the life of scanners and ensure smooth operation.

## 246.2 Fluid Precautions

If fluid and/or moisture penetrates the seal surrounding the piezo scan tube, the microscope will be severely damaged. The piezo scan tube is protected by the piezo guard and a silicone o-ring. Atop the silicone o-ring sits a thin seal of sealant glue.


**Figure 246-1** Cross-section of scanner. Note the location of silicone glue seal.

### *Document Revision History: Support Note 246*

Rev.	Date	Section(s)	Ref. DCR	Approval
→ Rev. A	05SEP97	Released.	0189	CCF


The ring itself is only a support for the glue, which acts as the actual seal. If the glue separates from the wall of the piezo tube or piezo guard, the seal has been compromised and must be repaired before it is used with fluids.

DI recommends that the seal be checked on a regular basis. To check the seal around the piezo scan tube, hold the scanner body underneath a low-power optical microscope. Focus on the o-ring and the silicone glue. Using a pointed instrument (like the tweezers supplied by DI), press GENTLY on the top of the o-ring.

**Figure 246-2** Close up of o-ring seal. Press gently on top of o-ring to check the integrity of the silicone glue seal.

If the sealant separates from the wall of the piezo scan tube or the piezo guard, the seal should be repaired before using the scanner with fluids.

**Figure 246-3** If the silicone glue seal separates from the wall of the piezo tube OR the scanner body, the seal must be repaired.

The silicone glue seal may be repaired using **Dow Corning #3140 RTV Silicone Glue**, or a good quality silicone glue thinned to paintable consistency with acetic

acid (vinegar). Note that the glue seal may only be repaired a limited number of times before the sealant layer becomes too thick and interferes with the accuracy of the piezo tube. In this case, the entire o-ring must be replaced and another layer of silicone glue applied. O-rings are available from Digital Instruments upon request.

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**CAUTION!** Due to occasional problems some users have experienced with moisture wicking up through fluid cell seals, it is recommended that you avoid prolonged immersion in fluids. After completion of your experiments, always remove the cell from the fluid, detach it from the scanner, and dry thoroughly prior to storage. Any moisture present on the end of the scanner must be dried IMMEDIATELY to prevent shorting the piezo.

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