

How to calibrate a LynxEye or LynxEye XE/XET

Primary optics: 1mm divergence slit and 2.5 deg. soller (hidden under plate between slit and tube for D2 Phaser)

Sample: NIST Corundum (Bruker corundum for D2)

Secondary optics: 2.5 degree soller

Open DIFFRAC.Maintenance and go to the Detector tab.

Change the calibration settings to match the following screenshots.

For D2 (the theoretical peak should be typed in as 35.129 and it may change to a different value automatically):

Calibration by Measurement			
Theoretical 2Theta Peak	[°]	35.12661151	<input checked="" type="checkbox"/> Theta=2Theta/2
Theoretical Theta Peak	[°]	17.563305...	Use Actual
Detector Opening	[°]	5.856312812	
Measurement Range	[°]	3.000	
Step Size	[°]	0.300	
Time per Step	[s]	10.000	
Sample Rotations			
Phi	Speed [1/min]	15.0	

For D8 Advance/Discover:

Theoretical 2Theta Peak	[°]	35.1490	<input checked="" type="checkbox"/> Theta=2Theta/2
Theoretical Theta Peak	[°]	17.5745	Use Actual
Detector Opening	[°]	3.297358830	
Measurement Range	[°]	1.0000	
Step Size	[°]	0.1000	
Time per Step	[s]	10.000	
Sample Rotations			
Phi	Speed [1/min]	15.0	

DIFFRAC.DETECTOR - User: Lab Manager - Application Type: Powder Diffraction - Instrument: MeasSrv(D2-004583)/D2-004583

File Edit View Detector Help

TOOLS DETECTOR CONFIGURATION DB MANAGEMENT LOG

Detector Lynxeye

Load and Save the Detector State
Resolution And Offset
Results

Resolution and Zero Offset Determination

Calibration by Theoretical Calculation

Number of Channels		192
Channel Size	[mm]	0.075
Detector Size	[mm]	14.4
Secondary Track Radius	[mm]	141
Detector Angle	[°]	5.831266462

Calculate Advanced Settings

Calibration by Measurement

Theoretical 2Theta Peak	[°]	35.14661151	<input checked="" type="checkbox"/> Theta=2Theta/2
Theoretical Theta Peak	[°]	17.573305...	Use Actual
Detector Opening	[°]	5.823381089	
Measurement Range	[°]	1.000	
Step Size	[°]	0.100	
Time per Step	[s]	10.000	

Sample Rotations


Phi	Speed [1/min]	0.0
-----	---------------	-----

Save Experiment


Found New Hardware
Your new hardware is installed and ready to use.

start | My Computer | 12:43 PM


Click on the Measure button.



The calibration measurement will take a couple of minutes.



After the calibration measurement is complete you will get a satisfactory notice with a Green or yellow dot next to it. Now click Apply to Instrument.



When the dialog box comes up, press Yes. You will have to then reconnect to the instrument via Measurement Server-> Instrument status.

