## General Safety Instructions

- **Introductory User Manual**  
**Original Instructions (English)**  
**Translation of the Original Instructions (German, French, Spanish)**

The reproduction, transmission or use of this document or its contents is not permitted without express written authority. Offenders will be liable for damages. All rights reserved.

We have checked the contents of this manual for agreement with the hardware and software described. Since deviations cannot be precluded entirely, we cannot guarantee full agreement. However, the data in this manual are reviewed regularly and any necessary corrections are included in subsequent editions. Suggestions for improvement are welcome.

All configurations and specifications are subject to change without notice.

Order no. DOC-M90-ZXX006 V1. Update: Jan 26, 2010.

© 2009 - 2010 Bruker AXS GmbH, Karlsruhe, Germany.

All trademarks and registered trademarks are the sole property of their respective owners.

Printed in the Federal Republic of Germany.

Bruker AXS GmbH  
Östliche Rheinbrückenstr. 49  
76187 Karlsruhe, Germany  
Tel. +49 (7 21) 5 95-28 88  
Fax +49 (7 21) 5 95-45 87  
info@bruker-axs.de  
www.bruker-axs.com

# General Safety Instructions

<b>Original Instructions .....</b>	<b>1-1</b>
<b>Introduction .....</b>	<b>1-1</b>
About This Manual .....	1-1
Disclaimer and Liability .....	1-2
Warning symbols.....	1-2
<b>General Safety Precautions .....</b>	<b>1-4</b>
Qualified Personnel.....	1-4
Correct and Intended Usage.....	1-4
When System Arrives .....	1-5
In Case of Damage .....	1-5
Room Planning and Installation .....	1-6
Operation .....	1-8
X-ray Safety System .....	1-8
Resetting Safety System after alarm appears .....	1-9
<b>Hazards .....</b>	<b>1-11</b>
Radiation Hazards .....	1-11
Electrical Hazards .....	1-13
Mechanical Hazards .....	1-14
Machinery Safety .....	1-16
Hazards Caused by LED Illumination .....	1-17
Thermal Hazards.....	1-18
Hazards Generated by Materials and Substances .....	1-19
Hazards Caused by Gases .....	1-22
Hazards Caused by Magnetic Fields .....	1-22

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## Introduction

### About This Manual

This manual contains general information and guidelines about the Bruker AXS instruments, which have to be observed to ensure personal safety, as well as to protect the product. These notices are highlighted in this manual by the warning symbols and are marked as follows according to the level of danger.

For simplicity the types of instrument when mentioned in the text are referred to as XRD (X-ray diffraction) and XRF (X-ray fluorescence) instruments. These categories include: wavelength and energy dispersive XRF instruments, powder, materials research, or single crystal XRD instruments (also known as SCD), and small angle X-ray scattering instruments.

#### **Danger**

The word 'Danger' indicates that death, severe personal injury, or substantial property damage can result if proper precautions are not taken.

**Warning**

The word 'Warning' indicates that personal injury or significant property damage can result if proper precautions are not taken.

**Caution**

The word 'Caution' indicates that minor personal injury or property damage can result if proper precautions are not taken.

**Note**

The word 'Note' draws your attention to particularly important information on the product, handling the product, or to a particular part of the product documentation.

## Disclaimer and Liability

Bruker AXS has checked the contents of this manual for agreement with the hardware, firmware, and software described. Since deviations cannot be excluded entirely, Bruker AXS cannot guarantee full agreement. However, the information in this manual is reviewed regularly and any necessary corrections will be included in subsequent editions. Suggestions for improvements are welcome.

If at any time there is a conflict between the safety information contained in this manual and any relevant national, state or local regulations, the local rules always take precedence.

## Warning symbols

The following indicators and symbols may be displayed on the system. Strictly obey all instructions and warning text printed on the labels which are attached to the various parts of the equipment.

**Radiation Danger!**

**Live Part! Risk of Electric Shock!**

**Caution! Read the Operating Instructions!**

**Protective Ground/Earth Terminal**

**Follow all valid national, state, and local regulations for disposal.**

**Danger of Injury! Danger of Crushing!**

**Hot Surfaces!**

**Danger of Toxic Substances!**

**General Biohazard Warning!**

## General Safety Precautions

### Qualified Personnel

All Bruker AXS X-ray systems or parts of them must only be setup and operated in accordance with this manual. Only suitable trained and qualified personnel shall be allowed to work on this equipment.

The installation, maintenance and repair of the system may be carried out only by personnel who are authorized and trained by Bruker AXS.

All repairs, adjustments and alignments performed on any components of the system (including the host computer) must be carried out in strict accordance with the established safety practices and standards of the country in which the equipment is installed.

### Correct and Intended Usage

The instrument and its components may only be used for the applications described in the technical description or in the Bruker AXS literature (e.g. XRD/XRF LabReports or XRD/XRF Application Reports). You can request those reports from your local Bruker AXS Service or Sales Representative. It may only be used in conjunction with devices or components from other manufacturers which have been approved or recommended by Bruker AXS, for example for automation solutions.

**Warning**

The instrument can only function correctly and safely if it is transported, stored, setup, and installed correctly, and if it is operated and maintained as recommended by Bruker AXS.

## When System Arrives

Upon arrival, check goods immediately for any visible damages. Even if there is only suspicion of damage, don't issue an entire receipt for the carrier. Carrier or other parties (e.g. shipping company, railway authority, port authority) must be held liable in writing (letter or reserve) for any damages or loss.

Contact a claim surveyor if the claim amount seems higher than EUR 3.000,00. The needed address is shown in the certificate of shipment or can be found under [www.vht-online.de](http://www.vht-online.de).

## In Case of Damage

In case of damage please send following information immediately, but no later than 5 days after receipt, to:

Bruker AXS GmbH

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76187 Karlsruhe

Germany

Fax: +49 721 595 4587

Phone: +49 721 595 2888

Email: [info@bruker-axs.de](mailto:info@bruker-axs.de)

subject: "transport damage"

- Supplier's reference no.
- Short description of the incident
- Estimated amount of the claim
- Photos
- Report of claim surveyor, if applicable

**Note**

Don't wait longer than 5 days with the announcement of damages even if you don't have the above mentioned information! Don't scrap, return or repair any damaged goods and/or packaging until further advice! In case of non-observance insurance coverage may expire.

## Room Planning and Installation

**Caution**

Carefully read this manual before installation and start-up of your Bruker AXS system. Also it is indispensable to read the 'Introductory User Manual' of the system as it contains valuable hints for a safe installation and operation of your instrument.

Careful room planning and the completed pre-installation are prerequisites for rapid final assembly and commissioning of your Bruker AXS system. Room planning and pre-installation has to be carried out with the assistance of the appropriate Bruker AXS office after suitable consultation.

The installation of your Bruker AXS system must be planned so that the instrument is easily accessible and that cooling air can flow without restrictions.

The floor or the platform (table for some systems) where your instrument is placed has to be leveled and needs to have a suitable load capacity.

Avoid placing the instrument in the direct solar exposure.

**Warning: Water Pressure**

For systems which need external water supply, the possible water pressure has to be restricted to a maximum value. If the local water pressure exceeds this maximum, then a water pressure reducer must be installed!

Usually the water inlet valves are only opened to the degree necessary. At times they might be entirely closed. As a result, water consumption is reduced to a minimum.

Assure that all tubing and all fittings outside the instrument are properly dimensioned and withstand the water pressure. Check the rating of the components before connecting the instrument.

For details on the maximum allowed pressure and further requirements on cooling fluids please refer to the 'Introductory User Manual' of your instrument.

**Warning: Gas Operation (P-10, Helium or Nitrogen)**

Systems with gas flow devices require external gas cylinders, if no in-house gas supply system is available. The cylinders must be secured so that the cylinders are prevented from falling over. Correct and secure installation of gas cylinders must be ensured at all times!

Follow applicable national, state, and local regulations for installation and operation of P-10 (argon-methane gas mixture), helium, and nitrogen gas operation and gas cylinders.

Provide a separate vent to open air to let off the P-10 gas!

For further information on special room planning requirements for your Bruker AXS system like wiring diagram for mains supply, connection of cooling water, technical data, temperature range, transport, storage and other details please refer to the 'Introductory User Manual' of your instrument.

## Operation

For any details about the operating elements, start-up and operation of your Bruker AXS system please refer to the corresponding 'Introductory User Manual' and 'User Manual' of your instrument.

## X-ray Safety System

All types of Bruker AXS instruments feature a failsafe X-ray safety system and are in compliance with the requirements for X-ray systems which are certified as a so-called "fully protected instrument" ("Vollschutzgerät") according to the German X-ray safety standards called "Röntgenverordnung".

Many instrument types also comply with the French radiation safety standard AFNOR NFC 74-100.

For details on valid certifications for other countries please refer to the 'Introductory User Manual' of your instrument or contact your local Bruker AXS representative.

The failsafe X-ray safety system of your Bruker AXS instrument basically consists of two independent and redundant safety lines. Both of them have to be faultless to enable high-voltage to switch on the generation of X-rays.

In case of any error on one of the safety lines, the X-ray radiation will be switched-off instantly.

To get a more detailed overview on the safety concept of your instrument please refer to 'Introductory User Manual'.

**Danger: X-ray Radiation**

Any manipulation on the safety circuits could cause severe injury and even lethal diseases due to X-rays.

Do not manipulate any interlocks or other safety-relevant devices! Any forced or irregular interlock manipulation will revoke the certification as a “fully protected system”.

Safety screws or screws with sealed heads serve as protection against unauthorized removal of X-ray safety related components. It is strictly forbidden for any person to remove parts which are protected by these types of screws. The positions of other safety related parts without safety screws or sealed screws are controlled by the safety system. Manipulations and modifications of the safety systems are strictly forbidden!

**Danger: X-ray Radiation**

If it seems that the safety system does not work correctly or that the radiation protection is impaired for any reason, you must switch off the system immediately and contact your local Bruker AXS Service representative. A frequent check of the proper operation of the safety system is recommended.

## Resetting Safety System after alarm appears

The X-ray Safety System is permanently checking the status of all components which are relevant for the X-ray safety. If there are any errors or inconsistencies the safety controller will shut down the X-ray generator immediately. The error which caused the shut-down will be captured and displayed in the appropriate software tools.

It is not possible – even after power down and restarting the instrument - to turn on the X-ray generator again until the safety system has been reset.

**Danger: X-ray Radiation**

The X-ray safety system may only be reset if the fault which caused the shut-down has been repaired and is no longer present! Contact your Bruker AXS Service for help.

The reset of the safety system is protected by password. For safety reasons the password itself must only be known by personnel being responsible for X-ray safety at the user's site. The reset procedure is described within the "Introductory User Manual" or the "User Manual" of your system.

The duties and responsibilities of the person in charge of the instrument are commonly described in local regulations and laws. Please contact your local authorities before starting to work with your X-ray analysis system.

# Hazards

## Radiation Hazards

### **Danger: X-ray Radiation**

XRD and XRF systems manufactured by Bruker AXS are analytical instruments with high-power X-ray sources. Shielding and safety equipment guarantee that the scattered radiation outside the instrument does not exceed the legal level during normal operation.

The direct beam of the X-ray source is extremely intense. Exposure to radiation for even a fraction of a second can cause severe burns. Longer exposure can cause severe injury. These injuries can lead to lethal diseases. Even without visible injuries or sensible pain, lethal diseases might evolve. This might happen even after many years. Strictly observe all national, state, and local regulations for the operation of X-ray systems!

Bruker AXS checks the emitted X-ray radiation of each single system prior to delivery under worst case conditions.

The instrument is in compliance with the requirements for X-ray systems which are certified as “fully protected instrument” (“Vollschutzgerät”) according to the German X-ray safety standards called “Röntgenverordnung”. The majority of the Bruker AXS instruments are also certified by the French radiation safety standard AFNOR NFC 74-100.

For the majority of instruments the enclosure of the XRF or XRD instrument serves as the X-ray radiation safety box. It must always be ensured that the enclosure and the setup of the goniometer correspond to one of the specified configurations and that the safety system is always in operation.

Manipulations and modifications of the safety systems are strictly forbidden!

### **Danger: X-ray Radiation**

All Bruker AXS X-ray systems are equipped with yellow or orange X-ray warning lights. Be careful! Whenever these are illuminated the X-ray tube may be powered up and generating X-rays.

**Danger: X-ray Radiation**

For the majority of instruments the enclosure of the XRF or XRD instrument serves as the X-ray radiation safety box. If any part of the enclosure (e.g. some systems feature front doors, which are made of a radiation protection glass) gets damaged by any reason you must immediately switch off the system. Contact your local Bruker AXS Service representative.

## Electrical Hazards

### **Danger: Live Parts! Risk of Electric Shock**

The system must be operated ONLY with such mains supply voltages as are listed on the instruments type plate.

As long as the system is connected to the external mains supply, some terminals of the internal components may present a danger of high-voltage shock.

In order to disconnect the system completely from the external mains network, the mains supply must be switched off at a wall socket (e.g. by removing the mains supply plug) or a switch external to the instrument.

### **Danger: Live Parts! Risk of Electric Shock**

Sometimes complex systems will be supplied by several power supply lines. The system and its components will be live until **all** power lines connected to the system are disconnected from the mains power supply.

### **Danger: High Voltage**

Inside the Bruker AXS XRD or XRF instruments voltages up to tens of thousands VDC may be present which however are not accessible from the outside. High voltages exist in the high-voltage generator and X-ray tube, as well as in the connecting high-voltage cable. In addition, the radiation detectors are operated with high voltages up to several kVDC.

If one of these high-voltage components is damaged, immediately disconnect the system from mains supply and contact your local Bruker AXS Service representative.

### **Danger: Electrical Hazards due to Leakage of Cooling Fluid**

Inside most of the Bruker AXS instruments are closed loop fluid circuits for cooling the X-ray tube, in some systems also for the X-ray generator and for the internal air-conditioning. In addition to that, some instruments are provided with water hoses which are connected to the laboratory's water supply system.

If there are any leakages (caused by faulty components or by accident during the refilling of the internal cooling fluid reservoir) the system must be turned off instantly! This must be done by switching off the automatic circuit breaker or power disconnecter which is located in the vicinity of the instrument. It is not sufficient to just turn off the system using the power-off button as some terminals inside the instrument will stay connected to mains supply voltage.

The system must stay powered-off until the leakage source has been repaired and all wet regions have been dried thoroughly. Be aware, there is a high danger of electrical shocks whenever there is leakage of cooling fluid!

## **Mechanical Hazards**

### **Warning: Moving Mechanical Components**

Some of the mechanical drives of the instrument are very powerful and can cause severe personal injury!

In order to prevent personal injury or damage, in standard operation all movements are stopped if the enclosure of the system is opened.

If you suspect that the safety system for motorized drives does not work correctly, immediately switch off the system and contact your local Bruker AXS Service Representative.

Manipulations and modifications of the enclosure control are strictly forbidden!

Several Bruker AXS instruments allow the user to change components like detectors, optics, sample handling devices etc.:

**Warning**

Before making substantial changes in mechanical setup, the power of the system must be turned off completely.

**Warning**

In order to prevent substantial damages in the system you must turn off the power of the system completely before connecting or disconnecting any cables attached to the X-ray tube, the radiation detectors, motors or to the various accessory components.

**Warning**

Do not place any objects which are not required for the operation of the system onto or into the instrument.

**Warning: Heavy Mechanical Parts**

Some components of the Bruker AXS systems are very heavy. There is a risk of injury when handling them during installation and maintenance.

Thoroughly observe the instruction procedures for installing the system and its accessories within the corresponding User's Manuals.

Heavy objects have to be lifted by two or more people together whenever possible. Use proper lifting techniques at all times. Strictly observe all national, state, and local regulations for permitted maximum weights. Also check the 'Introductory User Manual' of your instrument for details on heavy parts.

Use the following steps as an overview of proper lifting techniques.

1. **Plan:** practice the lift. While lifting, bend at the knees, keep your back straight, tighten your stomach and lift with your legs.
2. **Position:** keep your body close to the object you wish to lift—your stability increases the closer you are to the object. Keep your feet shoulder width apart.
3. **Movement:** avoid making awkward movements while holding a heavy object. Get help if the object is too heavy or cumbersome.

## Machinery Safety

Your Bruker AXS system fulfils the Machinery Directive 2006/42/EC of the European Union.

Accordingly, for all moving drives measures have been taken to avoid injuring people.

All drives used within Bruker AXS instruments have been categorized according to a so-called "Danger Class". When a drive cannot cause significant injuries it is considered as "safe" and is classified as class A. Dangerous drives are classified class B.

Class B drives can only be moved if the enclosure of the system is closed. Whenever the enclosure is opened, any class B drive will stop or not start to move.

In principle movements of class A drives are allowed even if the enclosure is opened. However within the Bruker AXS instruments in some applications class A drives are treated in the same manner as class B drives.

For details about the drives used in your system please refer to the 'Introductory User Manual' and 'User Manual' of your instrument.

**Warning**

Installation and servicing of drives must only be done by trained personnel. Do not connect motor cables to indexer board channels, which are not configured for the drive. This might damage the motor, the motor driver output stage or might lead to bad functionality of the drive. In rare cases motors might also lead to dangers if connected incorrectly.

## Hazards Caused by LED Illumination

**Warning: LED Light**

In some of the Bruker AXS systems the interior is illuminated by LED lights.

These are class 1 LED products. Do not stare into beam.

## Thermal Hazards

### Warning: Hot Surfaces

Several components (mainly all stepper motors or motorized components, the X-ray tube, optional cooling water pumps, vacuum pumps, air heating elements (XRF), optional heating stages ...) of the system get very hot during standard operation. They can cause severe burns when being touched. Under normal conditions these components cannot be accessed as they are covered by metal pieces. But, during maintenance and repair the hot surfaces might be touched by accident.

Be aware that the surface temperature of some internal components might exceed 72°C (172°F) whenever the system is switched on. Before touching any of these components turn off the instrument and wait until the components have cooled down sufficiently.

Don't change the configuration of the motors. Higher motor currents influence not only the motor's surface temperature but can also damage the motor and the control electronics.

### Warning: Cold Surfaces

Some components (e.g. optional low temperature stages of the XRD instruments) get very cold during standard operation. They can cause severe burns when being touched.

Observe the special operating instructions given in the manuals delivered with these accessories.

Before touching any of these components, turn off the instrument and wait until the components have heated up sufficiently.

## Hazards Generated by Materials and Substances

### Note: Interactions of Samples with X-rays

Before analyzing samples with your Bruker AXS instrument, make sure to choose the adequate measuring method. Before starting the measurement consider possible interactions of the sample with X-ray radiation.

### Danger: Analyzing Toxic Substances

To prevent any risk for health and safety the analysis of samples with XRD / XRF instruments, including handling and/or preparation steps, must be entrusted only to persons who have qualifications in accordance with their activities and functions.

Based on their technical training, their knowledge and experience and their knowledge of the applicable standards, these persons have to be able to evaluate the work assigned to them and recognize possible hazards. This also includes the knowledge of applicable accident prevention regulations, generally recognized safety regulations, EC guidelines and country-specific standards and regulations.

Carefully read and follow all safety and handling advises in the appropriate material safety data sheets.

When analyzing toxic substances, strictly follow national, state, and local safety regulations!

Only use sample holders which are recommended by Bruker AXS for the respective analyzing measurement.

### Warning: Analyzing Biological Active Substances

When analyzing biologically active substances, follow applicable national, state, and local safety regulations, along with IEC 61010-2-081.

**Warning: Beryllium**

Due to missing alternatives, a few components of an X-ray system require the use of parts made of Beryllium. Fumes or dust of Beryllium and its reaction products are compounds can be hazardous if inhaled or contacted! Corrosion of the Beryllium may occur during use of the instrument. Beryllium must not be cut, machined, or handled in any way.

Do not touch the X-ray tube's radiation outlet window, the detector's front window and (on some systems) the heating chamber's windows, as these windows could contain Beryllium. In case of an implosion of the X-ray tube or heating chamber's windows, you could get hurt by Beryllium fragments. This is why you must wear safety goggles and gloves when cleaning the instrument after such an incident.

Disposal of Beryllium must comply with all applicable national, state, and local laws and regulations.

**Warning: Disposal of instrument**

Follow all valid national, state, and local regulations for disposal and contact Bruker AXS Service for additional assistance.

The Instrument may contain Thallium (XRF, TLAP analyzer crystal), Beryllium (X-ray tube and detector) and batteries (some electronic boards).

**Warning: Battery**

Some of the electronic boards are equipped with batteries. Disposal of batteries must comply with all applicable national regulations.

**Caution: Oil**

Some systems (especially XRF instruments) require a vacuum pump. This vacuum pump uses oil as a lubricant. Disposal of oil must comply with all applicable national, state, and local regulations.

**Warning: Analyzing Harmful Volatile Substances**

When analyzing harmful volatile substances, strictly follow national, state and local safety regulations!

Only use the sample holders recommended by Bruker AXS for the respective measurements.

Provide a vent to open air for the exhaust gases of the vacuum pump!

**Warning: Thallium**

The analyzer crystal TLAP (XRF instruments) contains Thallium, an extremely toxic heavy metal. This crystal is water-soluble and can cause poisoning by skin contact.

Please handle this crystal with maximum care and avoid any skin contact! Strictly obey the applicable national, state, and local safety regulations.

Disposal of Thallium-containing parts must comply with all applicable national, state, and local regulations.

**Caution: Overpressure in the Instrument**

Fans provide the interior of your Bruker AXS system with fresh cooling air from the environment. On some instruments (especially the D2 Phaser, XRD) an overpressure is created inside which is relieved through the ventilation slots. When the enclosure is opened air flows from inside to outside of the housing. If the user intends to insert a powder sample it is possible that the powder is blown toward his face or body. In that case use only solid compressed powders and samples which can't be blown away from the sample holder.

## Hazards Caused by Gases

### **Warning: Use of Wrong Gas Type**

It is the responsibility of the user to assure that the correct gas type is connected to the instrument. Before opening the gas valves leading to the instrument, double check that the specified gas is used. In particular exclude the inadvertent use of any explosive or corrosive gas.

## Hazards Caused by Magnetic Fields

### **Caution: Magnetic Field**

Some of the Bruker AXS XRF instruments make use of strong magnetic fields (especially in the S2 RANGER). Please avoid placing discs, credit cards and any other magnetically susceptible items near the sample chamber.

## Hazards When Using a Temperature Chamber

### **Warning: Highly inflammable samples and cleaning agents**

Observe and adhere to your national safety regulations for handling the samples to be investigated (e.g. use of safety goggles, gloves, respiratory protection etc.).

Make sure that the setup location is sufficiently ventilated. The environment of the temperature chamber must be kept free of inflammable gases and vapours.

**Warning**

Samples can generate dust, fumes or gases during preparation, measurement and handling.

Opening of Temperature chambers can cause setting free dust, fumes and gases after heating and/or rapid flushing.

**Warning**

It is recommended to provide a fire extinguisher that is installed close to the analytical instrument.

## Maintenance of the System

To ensure that the system operates reliably, the maintenance intervals, as recommended by Bruker AXS, should be met.

Please refer to the 'Introductory User Manual' of your instrument or contact your local Bruker AXS representative for details on maintenance intervals and the required maintenance steps.

## Cleaning the System

### Warning

Before you start cleaning the instrument or components completely disconnect the system from mains supply!

For cleaning the side panels, the interior of the enclosure or the interior of the sample magazine, use dry cleaning utensils only. Do not use water or aggressive cleansing agents.

Some systems feature glass front doors that can be cleaned with any standard household-type window detergent. Use the same agent if your system comes with a transparent sample magazine cover.

For more detailed information on cleaning of your system please refer to the 'Introductory User Manual' of your instrument.

## Servicing Precautions

All repairs, adjustments and alignments performed at any component of the system (including host computer) need to be carried out by authorized and trained personnel in strict accordance with the approved work practices of the country where the equipment is installed only.

The following sections contain important information and safety hints especially for service personnel.

### Hazards Caused by Remote Network Access

Most Bruker AXS systems are equipped with an Ethernet network interface. Therefore, they can be controlled remotely by computers. Those can be part of the local company network or even a wide area network. The implemented safety features of Bruker AXS systems make sure that the instrument is – at all times – in a safe state during standard operation (with respect to X-rays, electrical voltages, and mechanical movements). In case of failure, generation of X-rays and movement of mechanical drives will be stopped immediately.

#### **Danger of X-ray Radiation, High Voltage and Moving of Mechanical Drives due to Remote Control via Network**

During maintenance or repair by service, it may be necessary to move motorized drives with enclosure opened, to perform X-ray leakage measurements or to bypass portions of the safety system.

In these situations, the instrument has to be physically disconnected from the network. This prevents unintentional remote control of the instrument via the network! In such cases, the instruments network interface must be connected directly to a computer located nearby.

This computer must also remain disconnected from any LAN networks during the entire duration of maintenance and repair work!

## X-ray Hazards

### **Danger: X-ray Radiation**

The direct beam of the internal X-ray source is extremely intense. Exposure to radiation for even a few seconds can cause severe burns. Longer exposure can cause severe injury. These injuries can lead to lethal diseases. Even without visible injuries or sensible pain, lethal diseases might evolve. This might happen even after many years. For this reason the radiation regulations must be strictly adhered to!

“Fully protected instruments” (“Vollschutzgeräte”) do not allow you to access the direct beam. However in servicing it might be necessary to work with an unprotected system while repairing a defective safety system or for alignment purposes.

Working with an unprotected system is dangerous. Alignment work with an unprotected unit must only be carried out by persons who are subject to personal dose measurements and medical surveillance. Only the lowest values for the high-voltage and tube current are allowed to be used. Use long alignment tools! Never hold your hand in the primary beam! Only remove those parts of the shielding whose removal is essential for the alignment.

**Never leave the instrument unattended in this mode!**

### **Note: On-site Check of the X-ray Safety System**

After any replacement of parts related to the X-ray safety system or any components of the X-ray generator it is necessary to verify the correct function of the X-ray safety system. Please refer to the ‘Introductory User Manual’ or the ‘Service Manual’ of your system for further details.

## Electrical Hazards

### **Danger: Accessible high voltages and mains supply voltages during servicing tasks**

Mains supply voltage (200-240VAC) will be present at the terminals and wires within the mains distribution unit. Additionally, many controllers (XRF or XRD instrument controller, detector controllers, etc.) are operated with mains supply voltage.

During maintenance and repair it might be required to operate the instrument while some panels are removed. In this case, you must be aware that many accessible terminals carry dangerous and possibly lethal supply voltages (for details see the mains wiring schematics). **Never leave the instrument unattended in this mode!**

Additionally, for error tracking it might be required to disconnect the X-ray tube's high voltage cable from the X-ray generator. In this case, please note that the high voltage output socket of the X-ray generator could carry very dangerous and possibly lethal high voltages. Strictly follow the safety instructions of the 'Service Manual'.

In addition to that, for servicing tasks, it could be necessary to disconnect the detector power supply cables from the high voltage outputs of the corresponding detector interface boards. The output currents of these high voltage supplies are internally monitored and limited to non-critical values. But nevertheless, do not touch any of these output sockets while the detector high voltage is turned on.

### **Danger: High Voltage**

After disconnecting the system from mains supply, high voltages are still present. They discharge over time. Therefore, wait at least 10 minutes before accessing any HV parts.

### **Danger: Live Parts! Risk of Electric Shock**

In some of the Bruker AXS systems the interior of the enclosure is illuminated by a standard fluorescent light tube located on the ceiling.

Turn off the system's power supply before removing/inserting the light tube.

## Hazards Caused by ESD (Electro Static Discharge) or Incorrect Mounting

The electronic modules contain electro statically sensitive devices. Before any components or modules are touched, the service person shall discharge himself / herself by touching an earthed object and take care that all ESD precautions are taken when working on the electronics.

The components must only be connected to the plugs of the modules provided for them. The block diagram (see 'Service Manual') provides information. The plugs of the modules are not protected against the insertion of incorrect components.

### Caution: Damage of Printed Circuit Boards

Before removing or plugging any cables from/to any printed circuit boards the system must be turned off by all means! Otherwise, the drivers' output stages and logical circuits can be destroyed severely.

### Danger: System Damage

Incorrect connection of modules or components may result in substantial damage in the system. Therefore thoroughly check all such operations before switching on the power supply.

### Note: Failure due to Incorrect Mounting of Printed Circuit Boards (only systems featuring a control rack)

It is very important that all printed circuit boards of the control rack are fixed correctly using **the designated screws**. Electrical noise and spikes coming from the cable shields must be conducted through the front panel to the metal rack of the controller.

If the screws are not fixed correctly then these noise signals will be conducted through the electrical circuits on the PCB and can cause strange behaviour of the system or sporadic functional problems.

## Hazards While Operating in Alignment or Servicing Modes

Some of the Bruker AXS instruments feature special alignment and servicing modes. These modes allow operation of mechanical drives while the enclosure of the instrument is open. Following safety precautions have to be thoroughly considered whenever the instrument is operated in these modes.

### **Danger: Injuries due to Moving Parts and Components**

Bruker AXS systems and their accessories contain a variety of motorized axes. During the normal operation mode these components are completely enclosed. But, during alignment and servicing there is an enhanced risk of injuries. For these modes the various drives might be accessible while moving. **Never leave the instrument unattended in this mode!**

The majority of the motors can be considered as safe drives that cannot cause any injuries at all when getting in touch due to their very limited power. But, there is a number of drives that can raise significant crushing, shearing, cutting or impact hazards to the operator while moving!

For more detailed information on the drives of your instrument please refer to the 'Introductory User Manual', 'User Manual' and 'Service Manual' of your instrument.

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