

## Electrochemical Cell for VeeMAX III™

The information in this publication is provided for reference only. All information contained in this publication is believed to be correct and complete. PIKE Technologies, Inc. shall not be liable for errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this material. All product specifications, as well as the information contained in this publication, are subject to change without notice. This publication may contain or reference information and products protected by copyrights or patents and does not convey any license under the patent rights of PIKE Technologies, Inc. nor the rights of others. PIKE Technologies, Inc. does not assume any liability arising out of any infringements of patents or other rights of third parties.

This document contains confidential or proprietary information of PIKE Technologies, Inc. Neither this document nor the information herein is to be reproduced, distributed, used or disclosed, either in whole or in part, except as specifically authorized by PIKE Technologies, Inc.

PIKE Technologies, Inc. makes no warranty of any kind with regard to this material including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Copyright 1991-2014 by PIKE Technologies, Inc., Madison, WI 53719. Printed in the United States of America. All world rights reserved. No part of this publication may be stored in a retrieval system, transmitted, or reproduced in any way, including but not limited to, photocopy, photograph, magnetic or other record, without the prior written permission of PIKE Technologies, Inc.

Address Comments to:

PIKE Technologies, Inc.  
6125 Cottonwood Drive  
Madison, WI 53719

Phone	(608) 274-2721
Fax	(608) 274-0103
E-mail	info@piketech.com
Web Site	www.piketech.com

Jan. 1, 2014

## Contents

Introduction	1
Cell Diagram	1
Unpacking Your Accessory	2
Packing List	2
Cell Assembly	3
Electrochemical Cell Installation	4
Cell Customization	5
Precautions	5
Replacement Parts and Options	5

## Introduction

PIKE Technologies' Electrochemical Cell is designed for in situ studies of redox processes and of electron transfer mechanisms. The Electrochemical Cell uses the VeeMAX III variable-angle specular reflectance/ATR VeeMAX platform. The cell is designed for ATR experiments. The angle of incidence may be varied to optimize spectral results. The ATR crystal is easily removed from the cell mount, which allows changing of the crystal type and facilitates the application of ATR surface coating. Additionally, the cell may be completely sealed for air-sensitive experiments. Electrodes are not included.





Figure 1. Cell Diagram

## Unpacking Your Accessory

In order for you to quickly verify receipt of your accessory, we have included a packing list. Please inspect the package carefully.


### Packing List

Electrochemical Cell Manual  
PN 350-001332  
Quantity 1



The image shows the cover of a manual titled "Electrochemical Cell for VeoMAX III". It features the PIKE logo and a photograph of the electrochemical cell.

Electrode Holder with Bracket  
Quantity 1



The image shows a black metal electrode holder with a bracket and a central electrode holder.

Teflon Bottom  
Quantity 1



The image shows a white, rectangular Teflon bottom plate with a central circular hole and four smaller holes at the corners.

Teflon Lid  
013-3201  
Quantity 1



The image shows a white, circular Teflon lid with four small holes around the perimeter.

Teflon Liquid Cup  
Quantity 1




The image shows a white, cylindrical Teflon liquid cup with a central hole and four small holes around the top edge.

Delrin Crystal Holder  
Quantity 1  
(Sold Separately)



The image shows a black Delrin crystal holder with a central vertical slot and four holes at the corners.

Screws  
4-40 x .25" (4) 10-32 x .25" (2)  
4-40 x .50" (4)




The image shows a collection of screws: four 4-40 x .25" screws, two 10-32 x .25" screws, and four 4-40 x .50" screws.

## Cell Assembly

1. With the liquid cup flipped over, center and place the ATR crystal over the small O-ring as shown below. The liquid cup is tapered. The larger inner diameter is the top of the cup and the smaller inner diameter is the bottom of the cup, where the crystal is placed.




2. Place the Teflon rectangular bottom over the ATR crystal and align with the screw holes on the bottom of the cup. The crystal orientation with the flat edge aligned perpendicular with the recess of the Teflon rectangular bottom. Using a 3/32" hex wrench, attach the Teflon bottom to the liquid cup using four 4-40 x .50" screws. Do not over-tighten the screws as this may strip the threads.



3. Place the Delrin crystal holder plate over the screws which were installed in step 2. The recess on the plate should be seated against the crystal. Secure using four 4-40 x .25" screws. Do not over-tighten the screws as this may strip the thread of the Teflon rectangular bottom. The Delrin crystal holder is customized for ATR crystal face angle. 45° and 60° Delrin holders are available.



*45° Delrin Holder*



*60° Delrin Holder*



4. It is recommended that integrity of the seal be checked prior to placing the electrochemical cell onto the VeeMAX III.

## Electrochemical Cell Installation

1. Included with the electrochemical cell is an electrode holder. If this is required for your experiment, attached the electrode bracket to the back of the VeeMAX III.



2. Remove the electrode bracket/micrometer by pulling upward below the micrometer to remove it from its seating pins.
3. Replace the specular insert located on the top of the VeeMAX with the assembled Electrochemical Cell.



## Cell Customization

PIKE's Electrochemical Cell for the VeeMAX III has been designed to be adaptable to custom experimental setups. Drill holes in the Electrochemical Cell Teflon lid to accommodate electrodes and purge lines. The micrometer holder located on the back of the VeeMAX III top plate may be used to set the distance between an electrode and crystal if required. The Derlin™ plug in the micrometer arm may be machined to accommodate electrode diameter.

Contact between an electrode and the crystal may be made through a user customized spring-type device. As an alternative, a thin foil may be threaded underneath the O-ring sealing the bottom of the liquid cup.

## Precautions

The O-rings shipped with the electrochemical cell are Viton. Be sure to check the compatibility of the O-ring to solvents being used in the electrochemical cell. Perfluoroelastomer O-rings are available. Contact PIKE Technologies to order.

Always check the integrity of the electrochemical cell seal prior to placing the electrochemistry onto the VeeMAX III base optics.

## Replacement Parts and Options

Part Number	Description
013-3110	ZnSe Crystal, 45°
013-3130	ZnSe Crystal, 60°
013-3112	Ge Crystal, 45°
013-3132	Ge Crystal, 60°
013-3114	Si Crystal, 45°
013-3134	Si Crystal, 60°
013-3115	ZnS Crystal, 45°
013-3135	ZnS Crystal, 60°
013-3345	45° Crystal holder
013-3360	60° Crystal holder





---

6125 Cottonwood Drive · Madison, WI 53719-5120 · (608) 274-2721 (TEL) · (608) 274-0103 (FAX)  
info@piketech.com · www.piketech.com